

# CRYSTALIZUAR DISH (W029)

Category: Lab Tools



**Product Code:** W029

**Category:** LAB TOOLS

**Sub-Category:** GLASSWARE

- 
- A crystallizing dish in engineering testing laboratories is used for evaporating solvents and facilitating the crystallization of dissolved substances. Typically made of glass or

porcelain, it allows for the gradual formation of crystals, often used in chemical analyses and purification processes.

---

## **ORDERING INFORMATION**

Item	Code
CRYSTALIZUAR DISH - Ø 140 mm	W029G140H

---

*Keywords: crystallizing dish, W029 crystallizing dish, ALFA Testing Equipment, laboratory glassware, solvent evaporation, crystal formation, engineering lab equipment*