

CRYSTALLIZING DISH (W029)

Category: Lab Tools



Product Code: W029

Category: LAB TOOLS

Sub-Category: GLASSWARE

-
- A crystallizing dish in engineering testing laboratories is used for evaporating solvents and facilitating the crystallization of dissolved substances. Typically made of glass or porcelain, it allows for the gradual formation of crystals, often used in chemical analyses and purification processes.

ORDERING INFORMATION

Item Name	Item Code
CRYSTALIZUAR DISH - Ø 140 mm	W029G140H

Keywords: crystallizing dish, W029 crystallizing dish, ALFA Testing Equipment, laboratory glassware, solvent evaporation, crystal formation, engineering lab equipment