

## MOTORIZED CEMENT FLOW TABLE (C-021/M)

- The test gives indications about the consistency of the cement mortars which is one of the most important characteristics for the mortars in the engineering aspects.
- The test is done by placing the specimen in the mould. After specifying the number of drops and running the machine, the plate will raise and drop along with the mould in a specific range for a specified number of times (entered by the user using the digital counter) depending on the test needs.
- Supplied with the testing mould.



Motorized Cement Flow Table (C-021/M)

### TECHNICAL SPECIFICATIONS

- Used for determining the consistency of cement mortars.
- Supplied with mould.