

CEMENT FLOW TABLE (C-021)

- The test gives indications about the consistency of the cement mortars which is one of the most important characteristics for the mortars in the engineering aspects.
- The test is done by placing the specimen in the mould. By turning the handle, the mould will raise and drop in a specific range for a specified number of times depending on the test needs.
- Supplied with the testing mould.



Cement Flow Table (C-021)

TECHNICAL SPECIFICATIONS

- Used for determining the consistency of cement mortars.
- Supplied with mould.

Code	Dimensions (± 1 cm)	Approximate Weight (kg)
C-021	34 x 26 x 32 (h)	13.8